



Maintenance of Cultures

Cultures sent to customers remain alive for at least 14 days provided that the following precautions are observed. Cultures should be unpacked immediately after receipt and stored at 15-18°C under low light intensity (north window, no direct sun light, or weak white fluorescent light). Screw caps or vessels should be loosened but not removed. Further maintenance or multiplication of cultures requires transfer into new culture media. This presupposes experience in simple microbial techniques.

Many species are cultivated and dispatched on agar media for safety reasons but develop their morphological characteristics only in liquid media, e.g. flagellates, colony-forming Volvocales and Chlorococcales. For teaching purposes these species should be transferred into liquid media 2-3 weeks before demonstration, e.g. into Soil Water Media, Basal Medium, or Desmidiacean Medium.

Culture Media

The following media have proved suitable for the maintenance of cultures in test tubes at the SAG for many years. The recipes originate from E. G. Pringsheim and W. Koch, unless stated otherwise. It must be emphasized that the maintenance medium indicated is not always the best medium for the cultivation of a species. There are other media which are just as suitable, e.g. those given in the catalogues of other culture collections of algae (Thompson *et al.*, 1988, Andersen *et al.*, 1991, Rippka & Herdman, 1992, Starr & Zeikus, 1993, Watanabe & Nozaki, 1994). Mass algal culture often requires more concentrated media (for recipes and methods consult Kuhl & Lorenzen, 1964; Starr, 1971; Stein, 1973; Guillard, 1975; Werner, 1982; Castenholz, 1988; Richmond, 2004; Andersen, 2005).

All solutions should be made up with de-ionized water. Media are usually prepared from stock solutions of macronutrients, trace metals, and vitamins which are added to a large proportion of the final volume of water in order to avoid precipitation.

Media may be used as liquid or solidified by 1.0-1.5% agar. Before sterilization the agar has to be dissolved in the medium in a steamer. After this test tubes should be filled with 10 ml of the hot medium, closed with cotton plugs, sterilized (usually by autoclaving at 121°C for 15 min.) and may be stored for several weeks, after cooling, in a refrigerator. Solid media for Cyanobacteria are prepared by mixing, after cooling to 50°C, equal volumes of separately autoclaved double strength solutions of the mineral salts medium and either agar to give a final agar concentration of 0.6-1.0 %.



References

- Andersen, R.A.; Jacobson, D.M. & Sexton, J.P. – Provasoli-Guillard Center for Culture of Marine Phytoplankton. Catalogue of Strains. 98pp. West Boothbay Harbor, Maine, USA, 1991.
- Andersen, R.A.; Morton, S. L. & Sexton, J.P. – Provasoli-Guillard National Center for Culture of Marine Phytoplankton. 1997 List of Strains. J. Phycol. 33 (suppl): 1-75, 1997.
- Castenholz, R.W. – Culturing methods for Cyanobacteria. In: L. Packer and A.N. Glazer, eds., Cyanobacteria. Methods of Enzymology 167 (1988), 68-93.
- Guillard, R.R.L. – Culture of Phytoplankton for feeding marine invertebrates. In: W.L. Smith and M. H. Chanley, eds., Culture of marine invertebrate animals. pp.29-60, Plenum Book Publ. Corp., New York, 1975.
- Kuhl, A. & Lorenzen, H. – Handling and culturing of *Chlorella*. In: D.M. Prescott, ed., Methods in cell physiology. Vol.1, pp. 152-187, Academic Press, New York and London, 1964.
- Rippka, R. & Herdman, M. – Pasteur Culture Collection of Cyanobacterial Strains in Axenic Culture. Vol.1, Catalogue of strains. 103pp., Institut Pasteur, Paris, France, 1992.
- Starr, R.C. – Algal Cultures – sources and methods of cultivation. In: A. San Pietro, ed., Photosynthesis. Part A, pp. 29-53, Methods in Enzymology vol. 23, Academic Press, New York, 1971.
- Starr, R.C. & Zeikus, J.A. – UTEX – The Culture Collection of Algae at the University of Texas at Austin. J. Phycol. Suppl. 29 (1993).
- Stein, J.R. ed. – Handbook of phycological methods. Culture Methods and growth measurements, pp. 448, Cambridge at the University Press, London, New York, 1973.
- Thompson, A.S.; Rhodes, J.C. & Pettman, I. – Culture Collection of Algae and Protozoa. Catalogue of strains. 164pp., Natural Environment Research and Council, England, 5th edit., 1988.
- Watanabe, M.M. & Nozaki, H. – NIES-Collection. List of strains, microalgae and protozoa. 4th edit., 127pp. The National Institute for Environmental Studies, Japan, 1994.
- Werner, D. – Biologische Versuchobjekte. Kultivierung und Wachstum ausgewählter Versuchsorganismen in definierten Medien. 432pp. Fischer Verlag, Stuttgart, New York, 1982.

Further recommended literature about culturing algae:

- Andersen, R.A. (ed.) (2005) Algal Culturing Techniques. Elsevier Academic Press, Burlington.
ISBN 0-12-088426-7.
- Belcher & Swale (1982) Culturing Algae - a guide for schools and colleges.
ISBN 1-871105-04-8 (ask for at ccap@sams.ac.uk). (Currently unavailable).
- Isaac & Jennings (1995) Microbial Culture. Bios Scientific Publ., Oxford.
ISBN 1-872748-92-9.
- Richmond (ed.) (2004) Handbook of Microalgal Culture. Biotechnology and Applied Phycology. Blackwell Publ., London. ISBN 0-632-05953-2.
- Streble & Krauter (2012) Das Leben im Wassertropfen (12. Aufl.). Kosmos (Franckh-Kosmos), Stuttgart, Germany. ISBN 13: 978-3440126349.
- Von Berg, Linne & Melkonian (2012) Der Kosmos-Algenführer (2. Aufl.). Kosmos (Franckh-Kosmos), Stuttgart, Germany. ISBN 13: 978-3440131732.

20. BG 11 Medium for Cyanobacteria (= BG 11)

Rippka, R. & Herdman, M. (1992), modified

	stock solution [g/100 ml]	nutrient solution [ml]
NaNO ₃	15	10
K ₂ HPO ₄ · 3H ₂ O	0.4	10
MgSO ₄ · 7H ₂ O	0.75	10
CaCl ₂ · 2H ₂ O	0.36	10
citric acid	0.06	10
ferric ammonium citrate	0.06	10
EDTA (dinatrium-salt)	0.01	10
Na ₂ CO ₃	0.2	10
micronutrient solution *		1
de-ionized or distilled water		919

* Composition of the micronutrient solution (from Kuhl and Lorenzen 1964):

Add to 1000 ml of de-ionized or distilled water:

H ₃ BO ₃	61.0 mg
MnSO ₄ · H ₂ O	169.0 mg
ZnSO ₄ · 7H ₂ O	287.0 mg
CuSO ₄ · 5 H ₂ O	2.5 mg
(NH ₄) ₆ Mo ₇ O ₂₄ · 4H ₂ O	12.5 mg

a) BG-11 Medium without Sodium Nitrate (=BG11-NaNO₃)

Prepare BG-11 Medium (Medium 20) without Sodium Nitrate (NaNO₃) and add 929 ml instead of 919 ml of water.

References:

Kuhl, A. & Lorenzen, H. (1964) Handling and culturing of *Chlorella*. In: D. M. Prescott, ed., Methods in cell physiology. Vol. I, p. 152-187, Academic Press, New York and London, 1964.

Rippka, R. & Herdman, H. (1992). Pasteur Culture Collection of Cyanobacterial Strains in Axenic Culture. Catalogue & Taxonomic Handbook. Vol. 1. Catalogue of Strains. 103 p., Institut Pasteur, Paris, France.

Stanier, R. Y., Kunisawa, R., Mandel, M. & Cohen-Bazire, G. (1971). Purification and properties of unicellular blue-green algae (order Chroococcales). Bacteriol. Rev. 35: 171-205.